

# OrthoLucent™ Hohmann Retractors

Carbon fiber material is lightweight, completely radiolucent, has strength equal to stainless steel, and can be steam sterilized

Also helps to prevent from scratching component surfaces.


#### PRODUCT NO'S:

4535-R [Modified]  
Blade Width: 18mm  
Overall Length: 10"

4550-R [Modified Blunt]  
Blade Width at Widest: 24.5mm  
Overall Length: 10.75"

4558-R [Standard]  
Blade Width: 16mm  
Overall Length: 9.625"

Ortho Lucent™ is a trademark of Ortho Tool, Inc.



Modified Hohmann

Modified Blunt Hohmann

Standard Hohmann

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362