

Ortho Mallets

with Easy Grip Handles

Designed to exert force on osteotomes, chisels, gauges and other instruments

These solid stainless steel mallets each have a comfortable 4½" grip made of a textured silicone that helps prevent the surgeon's gloved hand from slipping and helps maintain a solid grip.

Product No's:

7805
[Cloward-Style]
Overall Length: 8"
Handle Length: 4.5"
Head Weight: 1 lb.
Head Diameter: 1"
Ball Pean Diameter: .75"
Designed by J. Stiehl, MD

7810
[Small]
Overall Length: 8"
Handle Length: 4.5"
Head Weight: 1 lb.
Head Diameter: 1.3125"

7815
[Large]
Overall Length: 8"
Handle Length: 4.5"
Head Weight: 1.75 lb.
Head Diameter: 1.5"



Now with a larger diameter handle

textured soft silicone handle



Comfortable grip helps prevent the surgeon's gloved hand from slipping and helps maintain a solid grip.

FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

