

Ortho Mallets *with Easy Grip Handles*

Designed to exert force on osteotomes, chisels, gouges and other instruments

These solid stainless steel mallets each have a comfortable 4½" grip made of a textured silicone that helps prevent the surgeon's gloved hand from slipping and helps maintain a solid grip.



Easy Grip
Textured Soft Silicone
Handles

PRODUCT NO'S:	
7810 [Small]	
Overall Length: 8" (20,3 cm)	
Handle Length: 4.5" (11,4 cm)	
Head Weight: 1 lb. (.45 kg)	
Head Diameter: 1.3125" (3,3 cm)	
7815 [Large]	
Overall Length: 8" (20,3 cm)	
Handle Length: 4.5" (11,4 cm)	
Head Weight: 1.75 lb. (.8 kg)	
Head Diameter: 1.5" (3,8 cm)	



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362