

OrthoVise™ - Large


U.S. Patent #D398,208

See Reverse Side for
Small OrthoVise™ Information

- ▶ Made of stainless steel
- ▶ Models equipped with attachment bolts allow a slap hammer to be attached to the end, as well as to either side of the large OrthoVise™ (except Bent Jaw models), for greater adaptability
- ▶ Bent Jaw models are not available with side attachment bolts, but have an end attachment nut to accept a Standard Slap Hammer (#3925)
- ▶ A different size slap hammer is used for the large and small sizes of OrthoVise™
- ▶ Slap Hammers are designed with a hammer plate for the additional use of a mallet if desired

Side Attachment Bolts

Allows a Large OrthoVise™ Slap Hammer (#3950) to be attached to the side of the device. Available on Standard Large and Long Nose Large OrthoVise™ with Attachment Bolts only. (Cannot be added on later.)



STANDARD LARGE

PRODUCT NO'S:

OrthoVise™ Length: 10" (25,4 cm)

3980 with Attachment Bolts (two sides & end)
with Large OrthoVise™ Slap Hammer (#3950)

3980-01 with Attachment Bolts (two sides & end)
without Slap Hammer

3981 without Attachment Bolts
without Slap Hammer
with End Attachment Nut that accepts a
Standard Slap Hammer (#3925)



LONG NOSE LARGE

PRODUCT NO'S:

OrthoVise™ Length: 12" (30,5 cm)

3965 with Attachment Bolts (two sides & end)
with Large OrthoVise™ Slap Hammer (#3950)

3965-01 with Attachment Bolts (two sides & end)
without Slap Hammer



LONG NOSE LARGE BENT JAW

PRODUCT NO'S:

OrthoVise™ Length: 11.5" (29,2 cm)

3966 with Attachment Nut (end)
with Standard Slap Hammer (#3925)

3966-01 without Slap Hammer
with Attachment Nut (end) that accepts a
Standard Slap Hammer (#3925)



SLAP HAMMERS FOR LARGE ORTHOVISE™ MODELS

PRODUCT NO'S:

3950 [Slap Hammer for Large OrthoVise]
For use with 3965's, 3980's, 3981
Overall Length: 16.5" (41,9 cm)

3925 [Standard Slap Hammer w/16" Rod]
For use with: 3966's
Overall Length: 16" (40,7 cm)




THREADED ADAPTERS FOR LARGE ORTHOVISE™ MODELS

PRODUCT NO'S:

3980-02 [Small Adapter]
Changes Male End of a Slap Hammer to Female

3980-03 [Threaded Adapting Screw—Large]
For use with 3965's, 3966's, 3980's, 3981

Small Adapter   Female/Female Adapter Converts from Male/Male  Allows a Standard Slap Hammer (#3925) to be used with any Large OrthoVise™ with Attachment Bolts

Threaded Large  Can be used to append the Large OrthoVise™ with an Attachment Bolt for use with a Large OrthoVise™ Slap Hammer (#3950)

MADE EXCLUSIVELY
FOR INNOVATED BY
GERMANY

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362



OrthoVise™ - Small


U.S. Patent #D398,208

See Reverse Side for
Large OrthoVise™ Information

- ▶ Made of stainless steel
- ▶ Models equipped with attachment bolts allow a slap hammer to be attached to the end of the OrthoVise™ for greater adaptability
- ▶ Slap Hammer is designed with a hammer plate for the additional use of a mallet if desired
- ▶ A different size slap hammer is used for the large and small sizes of OrthoVise™

STANDARD SMALL

PRODUCT NO'S:	
	OrthoVise™ Length: 8" (20,3 cm)
3985	without Attachment Bolt without Slap Hammer
3985-01	with Attachment Bolt (end) with Small OrthoVise™ Slap Hammer (#3955)
3985-T	with Attachment Bolt (end) without Slap Hammer



LONG NOSE SMALL

PRODUCT NO'S:	
	OrthoVise™ Length: 9.5" (24,1 cm)
3975	without Attachment Bolt without Slap Hammer
3975-01	with Attachment Bolt (end) with Small OrthoVise™ Slap Hammer (#3955)
3975-T	with Attachment Bolt (end) without Slap Hammer



SLAP HAMMER FOR SMALL ORTHOVISE™ MODELS

PRODUCT NO:	
3955	[Slap Hammer for Small OrthoVise] For use with: 3975's, 3985's Overall Length: 8.75" (22,2 cm)



THREADED ADAPTER FOR SMALL ORTHOVISE™ MODELS

PRODUCT NO:	
3985-03	[Threaded Adapting Screw–Small] For use with: 3975's, 3985's



MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



www.innomed.net



INNOMED

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.