

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Stainless Steel OrthoVise with Slap Hammer

Designed with the option of using a slap hammer on three locations on the OrthoVise™ for better adaptability. The slap hammer can be attached to the end of the OrthoVise™ or on either side. The slap hammer is designed with a hammer plate if the additional use of a mallet is desired. Also available without the slap hammer in various sizes.

U.S. Patent #D398,208

Product No's:

Product No's:
Standard

- 3980 [Large (10") with Slap Hammer]
- 3981 [Large 10" without Slap Hammer or Attachments]
- 3985 [Small (8") without Slap Hammer or Attachments]
- 3985-01 [Small (8") with Slap Hammer]

Long Nose

- 3965 [Large (12") with Slap Hammer]
- 3966 [Large Bent Jaw w/Slap Hammer]
- 3975 [Small (9.5") without Slap Hammer or Attachments]
- 3975-01 [Small (9.5") with Slap Hammer]



**NEW STYLE
AVAILABLE**

*Bent Jaw
Long Nose
OrthoVise*

