

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2004 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net




FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION


Percussion Awl

Metal handle. 22cm total length.

Product No:
T1012



Used to begin tunneling the cortical bone during rotator cuff repairs or Bankart procedures



Crochet Hook Suture Passer

Product No:
T1013



Notched at tip to hook looped sutures and pull through tunneled bone

