

# Pin Inserter



*Used for 1/8" (3,2 mm) diameter pin insertion*

Designed to hold onto a 1/8" (3,2 mm) diameter pin while it is being inserted into a cutting block during total knee surgery or other applications where a 1/8" (3,2 mm) diameter pin is used. The pin inserter holds the pin tightly, yet releases it easily after insertion. It may be used with round or triangular end pins.

PRODUCT NO:

4020

Overall Length: 5" (12,7 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362