

Powers Femoral Tibial Space Gauge/Shim

Designed by Mark Powers, MD

Designed to help evaluate flexion/extension gaps, and as a shim for functional and kinematic alignment in total knee replacement



Profile

The device is steam sterilizable.

PRODUCT NO:

1193

Overall Length: 6.5" (16.5 cm)

Width: .75" (1.9 cm)

Thin End Depth: 1.5 mm

Thick End Depth: 2.5 mm



USA MADE

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



www.innomed.net

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362