

Rotator Cuff Awl

Designed to create a hole in the proximal humerus for reattachment of the rotator cuff

Ergonomic Radel® R Handle provides comfort and a secure grip.

Product No: S0107
Overall Length: 9.625"
Handle Length: 5"
Awl Prong Depth: 12.7mm



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

