

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2004 Innomed, Inc., All Rights Reserved.

**Innomed**

103 Estus Drive, Savannah, GA 31404  
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION

## Kodros Radiolucent Awl

Used with an image intensifier to locate holes in interlocking nails. Designed to help keep the surgeon's hands outside the image intensifier's rays, it is fully autoclavable. Two nail diameter sizes available.

Designed by  
S. Kodros, MD

Product No's: 8030-01 [3.7mm dia.]  
8030-02 [4.2mm dia.]

*Helps locate holes in interlocking nails*

