


# Radiopaque Goniometers

## Designed for Angle Determination

Transparent to X-ray—only white radiopaque markings show for easy reading. Used to check for X-ray distortion.

PRODUCT NO'S:
2000 [Set of 3]
2005 [Finger-size] Overall Length: 5"
2010 [Medium] Overall length: 8"
2015 [Large] Overall length: 14"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362