


# Retractor Clip for Smoke Evacuation Tube

Designed by James Saucedo, MD

*Repositionable stainless steel fastener designed to clip onto a retractor to help control the location of a smoke evacuation tube*

Allows for use on a 1/8" thick material with allowance for a "spring" fit.



**PRODUCT NO:**  
5466  
Dimensions: 27 mm Long  
18 mm Wide  
16 mm Deep



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

[www.innomed.net](http://www.innomed.net)



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

# INNOMED

**1.800.548.2362**