


Right Angled Subscapular Spreader – Blunt Tips

Designed by Edward McFarland, MD

Designed to hold the subscapularis muscle open when performing a subscapularis split approach to the glenoid



PRODUCT NO:

1652

Overall Length: 7.5" (19,1 cm)

Blade Depth: 2" (5,1 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362