

# Robb Cement Curette

Designed by William Robb, MD

*Designed to help remove cement around a knee or hip prosthesis*

Made of Delrin

**PRODUCT NO:**

5635  
Overall Length: 8"  
Freer End: 5mm  
Cup End: 10mm



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

**INNOMED**

1.800.548.2362