

Robb Cement Curette

Designed to help remove cement around a hip or knee prosthesis

Made of Delrin

Product No: 5635

Overall Length: 8"

Small Cup: 5mm • Large Cup: 10mm



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000

Fax 912 236 7766

Internet www.innomed.net

eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

