


Rogozinski Rongeur

Designed with cutting direction adjustments of 360°, allowing the instrument to be held in an ergonomic position for enhanced control, strength and precision


Designed by
Chaim Rogozinski, MD
Abe Rogozinski, MD



Includes an ejector pin for easy bone fragment removal. Rotating head, with no locks or preset positions eliminates the need to switch to an up-bite or down-bite rongeur.

70° CUTTING HEAD
OVERALL LENGTH: 12.5"
SHAFT LENGTH: 9.625"
Product No's:
5007-02 70°, 2mm Wide
5007-03 70°, 3mm Wide
5007-04 70°, 4mm Wide
5007-05 70°, 5mm Wide

90° CUTTING HEAD
OVERALL LENGTH: 12.5"
SHAFT LENGTH: 9.625"
Product No's:
5009-02 90°, 2mm Wide
5009-03 90°, 3mm Wide
5009-04 90°, 4mm Wide
5009-05 90°, 5mm Wide



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

