



Rogozinski Locking Needle Driver/Scissors

Designed by Chaim Rogozinski, MD

Designed with a quick lock & release handle, can drive a needle and cut a suture without changing instruments

PRODUCT NO'S:	
3083 [Standard]	
Overall Length: 6.5" (16,5 cm)	
3084 [Large]	
Overall Length: 7.75" (19,7 cm)	

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.