

Sandman Curved Bone Punch

Designed by Geoffrey A. Sandman, MD

Designed to help elevate a depressed tibial plateau fracture



PRODUCT NO:

5305

Overall Length: 14"

Shaft Length: 9.5"

Impactor Diameter: 12.5mm (.5")

MADE IN USA

FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362