



# Scott Uni & Total Knee Cement Removing Curette

Designed by Richard D. Scott, MD

*Sized, shaped and angled 90° to help with retrieval of posteriorly extruded cement behind the tibial component in both total and unicompartmental knee arthroplasty*

Ultra hard titanium nitride coating helps to extend curette life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

**PRODUCT NO:**

4247  
Overall Length: 9.625" (24,4cm)  
Overall Length: 5.25" (13,3cm)  
Cup Size: 4/0



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71



**1.800.548.2362**

© 2014 Innomed, Inc.