


Screw/Pin Removal Locking Pliers

Unique jaw designed to solidly grip and clamp onto a screw head, broken screw, or pin for removal



*Small
Jaw End & Bite*

*Standard
Jaw End & Bite*

Designed to securely grab pins as small as 1.4 mm (.055") up to 2.4 mm (.095")

Designed to securely grab larger pins, screw heads, or broken screws




Standard

Small

Reduced jaw size for smaller screws, pins and incisions

PRODUCT NO'S:	
S0142 [Standard]	Overall Length: 7.875" (20 cm) Jaw Width at End: 4 mm
S0142-01 [Small]	Overall Length: 7.875" (20 cm) Jaw Width at End: 4 mm



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362