

# Shark Tooth Grasper

Designed by Luis Ulloa

*Graspers sharp teeth help grasp onto tissue and bone*

Helpful in removing the labrum, and osteophytes around the acetabulum and around the glenoid. Also helps to remove meniscus, osteophytes and loose bodies. Helps facilitate working through a small incision without disrupting vision.



**PRODUCT NO:**

1798 [Standard]

Jaw Size: 6mm x 10mm

Overall Length: 10"

Shaft Length: 7"

1799 [Long Shaft]

Jaw Size: 6mm x 10mm

Overall Length: 12"

Shaft Length: 9"

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY



FREE TRIAL ON MOST INSTRUMENTS

© 2012 Innomed, Inc.

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362