

Short Tip Acetabular Retractor

Designed for retraction around the acetabulum

PRODUCT NO:

C1014

Overall Length: 13.25" (33,7 cm)

Blade Depth: 4.5" (11,4 cm)

Blade Width: 1" (2,54 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362