

Sidhu Tibia Clamp

Designed by Kuldeep Sidhu, MD

Designed to be used to securely grasp and remove an entire tibial wedge

The tapered lower pad slides under the cut tibial wedge without first having to use wedges, then, clamping allows the spikes in the upper pad to securely grasp the entire tibial wedge for easy removal.

PRODUCT NO:

3643

Overall Length: 10.25" (26 cm)

Pads: 60 mm x 30 mm

Spike Length: 7,5 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362