


# Sierra OrthoLucent™ EVA Pelvic Osteotomy Retractor

Designed by Rafael J. Sierra, MD

*Designed to help with retraction of the inner pelvis for direct visualization of the inner pelvis prior to iliac osteotomy, the retractor is made of a lightweight carbon fiber PEI composite material*



The OrthoLucent™ carbon fiber PEI material is radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.

**PRODUCT NO:**

4541

Overall Length: 10.7" (27,2 cm)

Blade Width at Widest: 24,5 mm

Blade Width at End: 10 mm

MADE EXCLUSIVELY  
FOR INNOMED IN  
SWITZERLAND



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**

# INNOMED