


# Soft Impact Mallets

*Now Available with Easy Grip Handles*

## **Provides shock-absorbing force**

Designed to have a shock-absorbing force, providing less bounce or wasted force. The Mallet is filled with a shock-absorbing media and has a flat striking surface to keep the Mallet centered on an instrument. The standard handle is manufactured of copolymer. The bottom can also be used to tap an implant in place.



### DIMENSIONS:

**2 lbs. Hammers**  
 Overall Length: 10.5"  
 Handle Length: 5"  
 Head Width: 3.5"  
 Head Diameter: 1.375"

**3 lbs. Hammers**  
 Overall Length: 11"  
 Handle Length: 5"  
 Head Width: 3.5"  
 Head Diameter: 1.875"

### Product No's:

7820 [2 lbs. with easy grip handle]  
 7832 [2 lbs. with delrin end and easy grip handle]  
 7837 [3 lbs. with easy grip handle]

**FREE TRIAL**  
 ON MOST INSTRUMENTS  
 CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

**Phone** 800 548 2362  
 912 236 0000

**Fax** 912 236 7766

**Internet** www.innomed.net

**eMail** info@innomed.net

**Innomed**  
 ORTHOPEDIC INSTRUMENTS  
 103 Estus Drive Savannah GA 31404

