

Bayonet Scissors


Long and slender for great viewing area

The tips are slightly curved for good visualization while cutting away tissue. The tip is angled 20° from the shaft and has a 7mm spherical radius cutout inside the tines.

PRODUCT NO:

H4020

Overall Length: 11.5"



Bayonet Disposable Surgical Blade Holder

Used mainly for MIS applications, the working area is extended to reach through most dilation tubes and retractors

For use with standard #11 & #15 disposable surgical blades.
Bayonet design allows for optimum visibility.

PRODUCT NO:

H4030



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362