


Bone Awl

Sharp narrow tip allows the surgeon a quick and easy solution to always placing a pedicle screw in the desired location regardless of the curvature of the spinus process

PRODUCT NO:
H5950



Bone Taps

Designed to provide a more secure channel for inserting pedicle screws



PRODUCT NO'S:
H5955 5.5mm
H5962 6.25mm
H5970 7.0mm
H5975 7.75mm
H5985 8.5mm



Disposable Drill

Cost effective 2mm diameter drill for cervical and trauma applications



PRODUCT NO:
H5945

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362