

# Disc Spreaders

Available in both straight and angled shafts, these spreaders are designed to properly size the vertebral disc space after decompression and before cage/bone insertion



## PRODUCT NO'S:

### **Straight Disc Spreaders**

H7107 7mm Straight

H7109 9mm Straight

H7111 11mm Straight

H7113 13mm Straight

### **Angled Disc Spreaders**

H7207 7mm Angled

H7209 9mm Angled

H7211 11mm Angled

H7213 13mm Angled

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362