

Distractor

Standard screw distractor is angled on the distal end and will work with most screw/rod systems up to 6.35mm diameter rods

PRODUCT NO'S:
H3710



Distractor with In-line Tips

This in-line version allows for distraction across the lamina thus eliminating any pressure or outside force applied to the pedicle screw construct

PRODUCT NO'S:
H3730



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED 

1.800.548.2362