



Plate & Rod Cutter

Designed to cut 3.5mm cervical rods and customized for various posterior cervical plating systems, all with one instrument

PRODUCT NO:
H3220



1/4" Rod Cutter with Depth Gauge

This all-purpose rod cutter has a patented guide within the jaws which helps minimize any twisting of larger rods, reducing possible breakage

When cutting in situ or non-in situ, the Depth Gauge, located on the jaw, indicates when the tips of the blade will clear the far side of the rod, ensuring the rod is completely cut through.

PRODUCT NO:
H3240



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362