

Rodholders

For use where excellent gripping of spinal rods are required without leaving any marks or possible fracture points. Available for any size rod. Distal ends may be straight or angled for better visibility or leverage

Rodholder – Angled Nose with Ratchet

PRODUCT NO:
H3500 3/16" (4.75mm) Angled Nose




Rodholder – Straight Nose with Ratchet

PRODUCT NO'S:
H3502 3/16" (4.75mm) Straight Nose
H3504 1/4" (6.35mm) Straight Nose




Rodholder – Angled Narrow Nose with Ratchet

PRODUCT NO'S:
H3506 7/32" (5.5mm) Narrow Nose
H3508 3/16" (4.75mm) Narrow Nose
H3510 1/4" (6.35mm) Narrow Nose



Rodholder – Straight Narrow Nose with Ratchet

PRODUCT NO:
H3512 7/32" (5.5mm) Narrow Nose with Ratchet



Rodholder – Angled Narrow Nose without Ratchet

PRODUCT NO:
H3514 7/32" (5.5mm) Narrow Nose without Ratchet



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED

1.800.548.2362