

Rod Pusher

PRODUCT NO:
H5990



Rod Persuader

Simple low cost solution for rod reduction with most pedicle screw systems

The tip is angled 20° from the shaft and has a 7mm spherical radius cutout inside the tines.



PRODUCT NO:
H5995
Overall Length: 11.5"

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362