




Screw Distractor

This unique instrument will distract a construct directly on the screw below the housing without effecting the angulation of the screw housing for proper alignment with the rod


PRODUCT NO'S:
H3720



Distractor with 180° Rotating Tips

Adjustable rotating tips allow for creating additional space for other instrumentation while distracting the disc space

PRODUCT NO'S:
H3740



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362