

Wide Cerebellar Retractors

General Gelpi-type retractors made of 304 stainless steel can be used for a variety of purposes

PRODUCT NO'S:	
H5300	30° Wide Cerebellar
H5310	60° Wide Cerebellar
H5320	90° Wide Cerebellar



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED 

1.800.548.2362