


Stanton Forward Retractors

Designed to work as "tissue pushers", helping to enhance exposure by allowing the surgeon or an assistant to push forward the opposite side of the wound while retracting (or not) the nearer side

Designed by
John L. Stanton, MD, FACS

Helps the assistant to hold the wound open without becoming fatigued and having to reach around to the opposite side of the wound to use another retractor.




Stanton Forward Ragnell Retractors



Shallow Forward Ragnell



Shallow Forward Ragnell



Product No's:
4510-01 [Shallow]
4510-02 [Deep]

Overall Length: 6.25"
Blade Offset from Handle: 1.625"
Blade Width: 3.1mm
Blade Depth: Shallow 13mm
Blade Depth: Deep 19mm


Stanton Forward Army-Navy Retractors



Shallow Forward Army-Navy



Deep Forward Army-Navy



Product No's:
4515-01 [Shallow]
4515-02 [Deep]

Overall Length: 9"
Blade Offset from Handle: 2.625"
Blade Width: 18mm
Blade Depth: Shallow 25mm
Blade Depth: Deep 43mm

Stanton Forward Senn Retractors



Shallow Forward Senn with Teeth



Deep Forward Senn without Teeth



Product No's:
4520 [Shallow With Teeth]
4525-01 [Shallow]
4525-02 [Deep]

Overall Length: 6.25"
Blade Offset from Handle: 1.625"
Blade Width: 9.5mm
Blade Depth: Shallow 13mm
Blade Depth: Deep 19mm

FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000

Fax 912 236 7766

Internet www.innomed.net

eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404