

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Subscapularis Spreader


Also used for retracting a split deltoid in mini-rotator cuff repairs.

15cm length.

Product No:
T1005



*Reaches deep to help
split the subscapularis
in a Jobe approach*



Kolbel Soft Tissue Retractors

Use facilitates the introduction of deeper retractors which are required for sufficient visibility of the glenoid, acromion and rotator cuff.

Product No's:
T1006 [Standard]
T1006-01 [Deep]



*Helps in the early phase to
retract soft tissue comprising
of the gleno-humeral joint*