



Bowen Suture Drills

***Used for suture placement during
arthroscopic Bankart repair***

Facilitates accurate penetration of Bankart soft tissue and placement of the suture through the scapular neck. Enhanced cutting surface proximal to the tip and recessed end facilitate suture delivery.

Product No's: (Sets of 2)

2300-01 [30 cm, 2.4 mm dia., 1 Slot, No Knurl]
2300-02 [30 cm, 2.4 mm dia., 1 Slot]
2300-03 [40 cm, 1.7 mm dia., 1 Slot]
2300-04 [40 cm, 2.4 mm dia., 1 Slot]
2300-05 [40 cm, 2.8 mm dia., 3 Slot]
2300-06 [36 cm, 2.8 mm dia., 4 Slot]
2300-07 [30 cm, 2.8 mm dia., 3 Slot]



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

