

TKR Knee Positioner

Provides stable positioning of the knee during surgery



Three (3) Sterile Pads/Wraps are included with each new purchase.

Leg Positioner Sterile Protective Pad & Wrap

Disposable, latex-free sterile foam pad and cohesive wrap helps protect patient from pressure sores, abrasions and possible neurological impairment while securing foot into the boot.

PRODUCT NO'S:	
2629-00	[Case of 10 Sets]
2629-L	[1 Set]



Each set includes 1 Pad and 1 Wrap



PRODUCT NO'S:	
2618-00	[Set]
Set Includes:	
2618-FP	[Aluminum Pivot Mount Footpiece] Length: 12.93" (32,8 cm) Height: 8.61" (21,9 cm) Width: 6.2" (15,7 cm)
2618-BP	[TKR Leg Positioner Base Plate] Length: 24" (61 cm) Height: .875" (2,22 cm) Width: 6.55" (16,6 cm)
2618-01	[TKR Foot Clamp Slide] Length: 3.75" (9,5 cm) Height: 4" (10,2 cm) Width: 4.375" (11,1 cm)
2629-L	[Sterile Pad/Wrap] Three included in set; product number is for one only
Optional & Replacement Parts:	
2629-00	[Case of 10 Sterile Pads/Wraps]



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.