

# Tendon Braiding Forceps

*Conically shaped fine point on one jaw separates tendon strands without causing long splits*

Very helpful when braiding tendons, the jaws are constructed to grasp and hold a small tendon stump securely. Features a delicate, sharp point and single pin tooth.

**PRODUCT NO:**

L5218 [Adult Size – Straight]  
Overall Length: 6" (150mm)

L5220 [Adult Size – Curved]  
Overall Length: 6" (150mm)

L1451-01 [Child Size – Straight]  
Overall Length: 4.75" (120mm)



**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



**INNOMED** 

**1.800.548.2362**