

Teurlings Medial Malleolar Clamp w/Wire Guide

Designed by Luc Teurlings, MD

Helps to stabilize the medial malleolar fragment during internal fixation



PRODUCT NO:

1803

Cannula Diameter: .062" (1.6mm)

Overall Length: 5.25"

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED 

1.800.548.2362