


Torx/Hex Adapter Set


Designed by Stephen M. Walsh, MD

Designed for conversion of a 3.5 mm screwdriver

Especially helpful when an articulated, universal joint driver is needed (i.e. acetabular screws)



Hex Bit to Torx Driver Adapter



Torx Bit to Hex Driver Adapter

PRODUCT NO'S:

8003-00 [Set – One Each]

Set Includes/Available Separately:

8003-01 [Torx Bit to Hex Driver Adapter]

Overall Length: .6" (1,54 cm)

8003-02 [Hex Bit to Torx Driver Adapter]

Overall Length: .6" (1,54 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362