

Universal Screwdriver Set

Helps eliminate the opening of multiple sterile packs when a specific size or style of screwdriver is needed

Helpful during revision total joint surgery where screws have been used, removal of bone plates, fracture fixation screws or bone graft screws. Set consists of: 6 (six) double ended screwdriver bits – small & large single slot, cross & cruciate, 3.5mm & 4.5mm hex, small & large phillips, small & large torx – a handle which accommodates any of the above bits, and a sterilization case.

NEW
BITS AVAILABLE



PRODUCT NO'S:

5195 [Complete Set with Storage Case] <i>Also sold individually</i>
5195-01 [Handle]
5195-02 [Straight (single slot)] <i>Large: 7x1.5mm, Small: 5x1mm</i>
5195-03 [Cross/Cruciate] <i>Large: 7mm, Small: 6mm</i>
5195-04 [Hex] <i>Large: 4.5mm, Small: 3.5mm</i>
5195-05 [Phillips] <i>Large: 4mm, Small: 3.5mm</i>
5195-06 [Small Torx: #10 & #15] New!
5195-07 [Large Torx: #20 & #25] New!



Set in Storage Case

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362