

Universal Deltoid Retractor

Designed for lateral retraction of the deltoid during arthroplasty or open shoulder procedures. Features two sizes for varying patient deltoid sizes.



PRODUCT NO:

1671

Overall Length: 8.9" (22,6 cm)
Large End Width: 2.75" (7 cm)
Large End Depth: 2.15" (5,5 cm)
Small End Width: 1.95" (5 cm)
Small End Depth: 1.75" (4,5 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362