

Vosburg Cannulated Periarticular Clamp

Designed by Caleb Vosburg, MD

Cannulated clamp tips allow passage of k-wires

By compressing the fracture with the clamp and then passing two k-wires, the clamp can then be removed to allow more working room and versatility when applying a plate.


PRODUCT NO:

1864

Overall Length: 13"

Handle Length: 8"

Ratcheted Opening from 2" to 3.5"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED 

1.800.548.2362