

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Weinraub Joint and Calcaneal Spreader

Designed by
Glenn M. Weinraub, DPM, FACFAS

Provides excellent joint exposure without blocking intra-articular or osteotomy access. Helps prevent slippage or falling out of the joint by placing the arms on either side of the area to be distracted, driving two pins and opening the joint.

Two models available:
one accepting up to
.063" (1/16" • 1.6mm) pins,
the other accepting up to
.110" (7/64" • 2.8mm) pins.

Product No's:
1870 [Up to 1/16" (.063") (1.6mm) Pin Diameter]
1872 [Up to 7/64" (.110") (2.8mm) Pin Diameter]



Designed to assist in the opening of small joints of the hand and foot for the application of fusion and graft techniques

