

# Wertz Anterior THA Femoral Elevator

Designed by Michael P. Wertz, MD

*Helps deliver the femur out of the incision during anterior total hip arthroplasty*

Inserted into the femoral canal for elevation, the knurled underside helps to reduce the chance of slippage.

**PRODUCT NO:**

6148

Overall Length: 12.625"

Blade Length: 3"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74



1.800.548.2362