

Whelan Femoral Neck Elevator

Designed by Edward J. Whelan, III, MD

Elevator has long tines to rest on the stronger bone at the base of the neck and calcar, and also fits well over the lesser trochanter and iliopsoas tendon for femoral broaching



PRODUCT NO:

3414

Overall Length: 13.75" (34,9 cm)

Depth from Bend: 1.5" (3,8 cm)

Blade Width: 2,4 cm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.