

# Wiater Shoulder Bone Hook


Designed by J. Michael Wiater, MD, FAOS, FAOA

*Large bone hook designed to retract the proximal humerus posteriorly to help provide exposure for glenoid reaming during open shoulder procedures*

Also useful for other large joint procedures.



**PRODUCT NO:**  
**5079**  
Curve Diameter: 2.68" (68 mm)  
Overall Length: 10" (25,4 cm)  
Handle Width: 4" (10,2 cm)



[www.innomed.net](http://www.innomed.net)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**