

Wilson Trigger Finger Retractor

Designed by Ralph V. Wilson, MD

Helps provide improved exposure, which frees up the assistant and gives easier and better visualization of the tendon sheath in a trigger finger procedure

Also useful in other small incision surgeries. Can use with Morton's neuroma and hand arthroplasty procedures.



PRODUCT NO:

1884

Overall Length: 4.25"

Blades: 6.5mm Wide x 12mm Deep

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED The Innomed logo consists of the word "INNOMED" in a bold, sans-serif font, followed by a stylized lowercase "i" that incorporates an American flag design.

1.800.548.2362