

**NEW**




# Wurapa Small Joint Compressor and Distractor

Designed by Raymond K. Wurapa, MD

*Designed to allow one-handed manipulation and deployment once fixation pins are placed*



Now available with two hole sizes on each instrument!



Designed to simplify several small joint procedures:

- ▶ Preparation of small bone non-unions before bone grafting and fixation
- ▶ Preparation of small joints for arthrodesis (e.g. partial wrist fusion)
- ▶ Useful for intercarpal stabilization while performing ligament reconstructions (e.g. scapholunate ligament repair/reconstruction)
- ▶ Distract and better evaluate small joints before determining final management.

For use with .045" (1.1mm) k-wire pins, which should be cut short above the pin guides to allow full access to the operative site.

PRODUCT NO'S:	
<b>Double 1.1mm &amp; 1.6mm Holes</b>	
1751 [Compressor] Compresses From: 28mm Overall Length: 4.625"	1752 [Distractor] Distracts to: 46mm Overall Length: 4.625"
<b>Single 1.1mm Hole</b>	
1753 [Compressor] Compresses From: 28mm Overall Length: 4.5"	1754 [Distractor] Distracts to: 46mm Overall Length: 4.5"

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362